








Accu-Chek® Guide Platform Meters and Test Strips

		
Meter	Accu-Chek Guide Me <i>Newest Accu-Chek meter</i>	Accu-Chek Guide
Weight <i>(incl. batteries)</i>	1.5 oz (43 g)	1.7 oz (48 g)
Size <i>(L x W x H)</i>	3 x 1.9 x 0.63 in (76 x 48 x 16 mm)	3.2 x 1.9 x 0.8 in (80 x 47 x 20 mm)
Display	LCD	Backlit LCD
Strip Ejector		●
Strip Port Light		●
Data Port	Bluetooth®, Micro USB	Bluetooth®, Micro USB
Languages	1	36
Lancing Device	Accu-Chek Softclix	Accu-Chek Softclix
Test Strip	Accu-Chek Guide	Accu-Chek Guide
Test Strip Size <i>(L x W)</i>	5 mm x 28 mm	5 mm x 28 mm
Application Area	 5mm <i>Widest application area among leading brands¹</i>	 5mm <i>Widest application area among leading brands¹</i>
Memory Test Results	720	720
Testing Time	< 4 seconds	< 4 seconds
Test Averages	7, 14, 30, 90 days	7, 14, 30, 90 days

¹Report of Dosing Area Measurement. Data on file. PR00017209.

The Choice is Yours.

Find the Accu-Chek® meter that fits best in your life.

		
<p>It is important for me to have the following features to assist me with the management of my diabetes:</p>	<p>Accu-Chek Guide Me <i>Newest Accu-Chek meter</i></p>	<p>Accu-Chek Guide</p>
<p>Access to Health-Related Apps</p>	<p>●</p>	<p>●</p>
<p>Meal Markers</p>	<p>Via mySugr app</p>	<p>On meter and mySugr app</p>
<p>Target Range Goals</p>	<p>Via mySugr app</p>	<p>On meter and mySugr app</p>
<p>Pattern Detection</p>	<p>Via mySugr app</p>	<p>On meter and mySugr app</p>
<p>Test Reminders</p>	<p>Via mySugr app¹</p>	<p>On meter and mySugr app¹</p>
<p>Comments</p>	<p>Via mySugr app</p>	<p>On meter and mySugr app</p>
<p>Features to Assist with Dexterity Challenges or Impaired Vision</p>	<ul style="list-style-type: none"> ▪ Large, easy-to-read display screen ▪ Larger dosing area² ▪ Easy-to-handle strip vial ▪ Proven ease of use³ 	<ul style="list-style-type: none"> ▪ Backlit display ▪ Test strip port light ▪ Strip ejector ▪ Larger dosing area² ▪ Easy-to-handle strip vial ▪ Proven ease of use³

¹Available with mySugr PRO

²Report of Dosing Area Measurement. Data on file. PR00017209.

³Usability evaluation of a blood glucose monitoring system with a spill-resistant vial, easier strip handling, and connectivity to a mobile app. Data on file. PDR_0000044114.

