

# The TRUE METRIX® Self-Monitoring Blood Glucose Systems are designed to help people confidently and accurately manage their diabetes.




## Select the Meter That Best Meets Your Individual Testing Needs

### TRUE METRIX® AIR – *Technology with Reach*

#### FEATURES:


- Bluetooth® Smart Connectivity
- No Coding
- Results as fast as 4 seconds
- 7-, 14-, 30-, 60- and 90- day averaging
- 6 Event Tags
- Stores 1,000 results with time and date



### TRUE METRIX® – *Accuracy Is Confidence*

#### FEATURES:

- No Coding
- Results as fast as 4 seconds
- 7-, 14- and 30-day averaging
- 6 Event Tags
- Stores 500 results with time and date



#### Featuring



Both systems use TRUE METRIX®  
Self-Monitoring Test Strips



Test Strips Made In The USA

Please visit [www.trividiahealth.com](http://www.trividiahealth.com) to find training videos on how to operate your meter, educational brochures, quick reference guides, and more. Much of our information is available in multiple languages.

**TRUE  
METRIX®  
AIR**  
SELF MONITORING BLOOD GLUCOSE SYSTEM

**TRUE  
METRIX®**  
SELF MONITORING BLOOD GLUCOSE SYSTEM

**TRIVIDIA**  
HEALTH®

NICO-4779 07/25 © 2025 Trividia Health, Inc. TRUE METRIX, TRIPLE SENSE TECHNOLOGY and the Trividia Health logo are trademarks of Trividia Health, Inc. TRUE METRIX AIR and TRUE METRIX are intended for self-monitoring blood glucose only and not for multiple patient use. Only TRUE METRIX PRO is intended for multiple patient use. All other product names and brand names are trademarks of their respective owners. All rights reserved.

# TRUE METRIX® Systems with advanced performance testing features to help you manage your diabetes.

Product Features	TRUE METRIX® AIR	TRUE METRIX®
Test Speed (Sec.)	As fast as 4	As fast as 4
Averaging	7/14/30/60/90	7/14/30
Ketone Test Reminder	Yes	Yes
Event Tags	4	4
Meal Tags	2	2
Sample Size (µL)	0.5	0.5
Alternate Site Testing	Yes	Yes
Download and Wireless Capabilities	TRUEmanager® Diabetes Management Software and TRUEMANAGER™ AIR app via Bluetooth® Smart	TRUEmanager® Diabetes Management Software
Testing Reminder Alarms (Sets)	Yes (4)	Yes (4)
Memory	1000	500
BG Range (mg/dL)	20-600	20-600
Hematocrit Range (%)	20-70	20-70
Blood Sampling	Fresh capillary whole blood from the fingertips or forearm	Fresh capillary whole blood from the fingertips or forearm

Visit



[www.trividiahealth.com](http://www.trividiahealth.com)

to learn more about the TRUE METRIX® AIR and TRUE METRIX® Systems!